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CONTINENTAL CONAIR LIMITED  
CHINA CREATIVE DEPARTMENT

Date	12/31/2021
Customer	VS-HK
Model No.	CWJ182H
Description	IB
Job No.	
Reference No.	IB-21/363
Revision No.	0
Color	1C+1C(Black)
Open Size	736mm(W) X 214mm(H)
Folded Size	92mm(W) X 214mm(H)
Material	80gsm Art/paper
Prepared by	Betty Bai
Approved by	Phoebe Laung/Andy Huang
Designer	
Design Job No.	
Vendor	JBY

折疊方法  
Folding method:  
16 pages (8 x 2)

CWJ182H

IB-21/363

## CONAIR PORTABLE USB RECHARGEABLE WATER JET

### INTRODUCTION

This portable water jet can improve your oral hygiene. It will effectively flush out the debris, reduce plaque between teeth, strengthen the gum and prevent tartar build-up.

### READ ALL INSTRUCTIONS BEFORE USING KEEP AWAY FROM WATER

### SAFETY INSTRUCTION: IMPORTANT SAFEGUARDS

When using an electrical appliance, especially when children are present, the following safety precautions should always be taken: Polythene bags over product or package may be dangerous. To avoid danger of suffocation, keep the wrapper away from babies and children. The bag is not a toy.

### SAFETY PRECAUTIONS

**WARNING:** Do not use or leave the power adaptor in a position where it can become wet or there is a danger of it falling into a basin or other vessels containing water. Keep the appliance dry, unless it otherwise specified.

This product is washable under water tap (sealing charging socket before cleaning).

The above mark indicates that this product is suitable for use in a bath or shower.


- This appliance is not intended for use by children under 5.
- Do not use any attachments other than those provided with this product.
- Use this appliance only for its intended use as described within this instruction booklet.
- Unless otherwise specified by a dental professional, fill the water tank with lukewarm water only. Do not overfill the water tank. It is recommended for the temperature of the water not to exceed 40°C to avoid scalding.
- Do not direct the water jet under the tongue, into the ear, or other delicate areas. The product can produce pressures that may cause damage in these delicate areas.
- Do not use the appliance without a water jet tip.
- Some mouthwashes may corrupt this product. Please do not add mouthwash to this product.
- Always unplug the appliance and seal charging socket when fully charged.
- This appliance is a personal care product. It is not suitable for commercial use or shared use by multiple patients.
- Do not drop or insert any object into an opening.
- Do not place on any surface while it is operating.
- Do not leave unattended while it is operating.
- Keep the cord away from heated surfaces.

- Never twist or pull on the power cord.
- After use, do not wrap the cord around the appliance as in time this may cause the cord to fracture. Coil cord loosely by the side of the appliance in storage.
- In the event of damage, discontinue use immediately and contact your dealer for instructions on returning it for examination or exchange.
- The charging voltage of this appliance is 5V DC. Before charging the appliance, make sure the output voltage of the charger corresponds to the voltage marked on the rating plate of the appliance.
- This appliance must be powered by safety extra-low voltage.
- Before cleaning, always unplug the appliance and secure the charging receptacle lid firmly.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure they do not play with the appliance.
- Do not use dishwasher or microwave oven to clean the jet tips, handle and USB cord.

### MEDICAL WARNING

- Do not use this water jet if you are fitted with a pacemaker or other implantable devices.
- If you have any medical concerns, consult your doctor before using this product.
- Patients with cardiac disease should consult their doctor before using this water jet.
- Contact your dentist prior to use if you have periodontal disease or have had oral or periodontal surgery or major orthodontic treatment in the previous 12 months.
- When you start any new oral hygiene program, your gums may initially bleed slightly. Contact your dentist if excessive bleeding occurs or if bleeding continues after two weeks of use.
- Do not use mouthwashes or oral rinses that contains iodine, decolorizer or tea tree oil.
- Remove any oral jewelry prior to use of this product.
- Do not use if you have an open wound on your tongue or in your mouth.
- Make sure only the tip of the water jet comes contact to the teeth, gum or tongue. If any other part of the handle contacts them, it may cause damage to these areas. Be careful not to press the tip against the gum, especially dental patients.
- This product is designed for oral cleaning only. Do not use the unit for any other purpose.

### PRODUCT DESCRIPTION


Application Scenarios of Water jet Tips		
Classic tip	Pressure water cleans food debris effectively	
Pik Pocket Tip	Ideal for removing debris in tooth pocket and prevent calculus	
Plaque Seeker	Ideal for retainer/implants/crowns/bridges cleaning	

### OPERATING INSTRUCTIONS CHARGING

**Note: The charging voltage of this appliance is 5V DC. Before charging the appliance, make sure the output voltage of the charger corresponds to the voltage marked on the rating plate (5V DC). A higher voltage will cause damage to the appliance.**

- If the charging indicator starts to flash when using or after use, the battery is low and needs recharging.
- Use the provided USB charging cable to connect the appliance to a computer USB port or a smart phone charger to charge the appliance (the charging voltage is 5V DC).
- Open the USB sealing plug on the bottom cover as indicated by the figure below, then insert the USB cable to charge.


- Once charging is started, the charging indicator on the front side of the handle will flash.
- The charging indicator will stay on when the appliance is fully charged.
- It takes 8 hours for the flusher to be fully charged (may vary depending on the environment).

### BEFORE USE

- Adjust the tip, then push the tip downwards until it is firmly attached to the handle (you will hear a click) (see the figure below). If you need to change the tip in the course of using, you should turn it off first. Simply press the changing button to change the tip.

**TIPS:** For best results, please replace the tip every 6 months. If the tip becomes loose or not "snapping into" the handle, you also need to replace the tip.


- Fill the water tank by either of the following methods:

-- Open the water tank lid, hold the water jet horizontally and fill with water. Replace the lid.

-- Remove the water tank, hold the water tank vertically and fill with water. Then replace the water tank from the vertical direction onto the main body (you will hear a click).


**NOTE:** Some mouthwashes may corrupt this product. Do not add mouthwash to this product.

- Press the power button to turn the appliance on. The LED indicator will light on.
  - Press the mode button to choose your desired mode.
    - There are 3 modes to choose from:
- Normal: suitable for normal cleaning of teeth.  
Soft: Suitable for users with sensitive teeth. Progressive power allows your gum to adapt to the pressure of the waterflow gradually, providing gentle cleaning without stimulation.  
Pulse: suitable for massaging the gum. This function can improve blood circulation and strengthen the teeth.
- Choose your desired mode by pressing the mode button before turning the power on.
  - This water jet can automatically memorize the previous mode selected by the user and prepare for the next use.

**TIPS:** In order to get a comfortable experience, when using the product for the first time, it is recommended to use soft mode for one week and change to normal mode until you get adapted to this product.

### USING THE PRODUCT


- Ensure the water jet is in off position. If it is operating, press the power button to turn it off.
- Press the mode button to select your desired mode. This water jet will automatically memorize the previous mode selected by the user. If your desired mode is consistent with the previous one, skip this step.
- Hold the water jet vertically, place the tip into your mouth with the tip directing at spaces between your teeth or the gum. Keep your mouth slightly close to avoid water spraying outside and water easy running out of mouth.

- Press the power button and hold for about 0.5 second to turn on the water jet to flush.

### NOTE:

- Direct the water jet at a right angle to the teeth or gum.
- Move along the teeth, clean the teeth in sequence.
- Direct the tip to the gumline at a right angle to the gum.
- Clean the teeth, braces, crowns and dentures thoroughly.
- Do not direct the water jet into tooth pocket.

- To turn the water jet off, simply press the power button.
- Automatic Power Off: This product is equipped with a 2-minute timer. When the appliance has not been operated for 2 minutes, it will automatically power off. To continue to


use, simply press the power button.

- For the sake of your health, please follow the following steps after each use:
  - Open the water tank lid and pour out the water from water tank.
  - Press the power button. The water jet will start to work and drain out the water in the appliance.
  - Press the power button again. The flusher will shut off.
  - Replace the water tank lid.
  - Wipe off any water on the flusher with a clean and dry cloth.

### DAILY CLEANING

#### NOTE:

- Before cleaning, always unplug the appliance and sealing charging socket with charging cap on the bottom of the main body. If charging cap is not covered properly, the appliance may be unable to charge or work after cleaning.**

- Use clean water or neutral detergent to clean this product. Do not use corrosive or abrasive cleansers because they would cause damage to the product. Do not clean this product with hot water over 50°C.
- Do not use dishwasher or microwave oven to clean the tips, handle and USB cable.
- If the residues in the tip or water tank are not removed properly, there may be a hygiene problem.

### CLEANING MAIN BODY:

- Clean the main body with a dry or damp cloth.
- Never immerse the main body into water.

### CLEANING WATER TANK:

- Wipe the water tank with a dry or damp cloth, or remove it to rinse under running water.
- If you are not going to use it for an extended period, please pour out the water and let it dry naturally.

### CLEANING WATER JET TIPS:

- Wipe with a dry or damp cloth, or remove it to rinse under running water. Then let it dry naturally.
- Do not bend, pull or twist the water hose.
- For best results, please replace the tip every 6 months. If the tip becomes loose or not "snapping into" the handle, you also need to replace the tip.

### STORAGE

If you are not going to use the appliance for an extended period, in order to avoid bacteria breeding, please follow the following steps:

- Open the water tank lid, pour out the water.
- Press the power button, the water jet will start to work and drain out the water in the appliance.
- Press the power button again to turn off the flusher.
- Remove the water tank by sliding it down.

- Clean the handle, water tank, tips with a clean dry or damp cloth.
- Leave the handle, water tank, tips and charging base dry out of direct sunlight.
- Close the water tank lid and slide the water tank back to handle.
- Store the handle, tips and USB charging cable in a cool dry place.


### DISPOSAL

This product is driven by rechargeable battery. Please remove the battery before discarding your product at the end of its life; and recycle or dispose the battery according to your local waste regulations. If you have any questions, please consult your local waste disposal office.

Note: Unless you are discarding this product, please do not attempt to disassemble it. This will cause damage to the product.

- Before removing the battery, ensure the battery has run out. If you need to run the battery out, press the power button and hold for 0.5 second to turn on. Run the flusher until the it stops operating. Repeat several times until the battery runs out completely.
- Before disassembling the unit and removing the battery, make sure the power is disconnected.
- The batteries must be disposed safely.
- Never replace the battery with other types of battery, especial non-rechargeable batteries. This could cause fire or an electric shock.
- Never throw batteries away with the normal household waste, make sure batteries to be disposed at an official recycling center or collection point.

### REMOVE THE BATTERY ACCORDING TO THE FOLLOWING STEPS:

### TROUBLESHOOTING

Fault Diagnosis		
Problems	Possible Reasons	Solutions
The appliance does not work.	The water jet has not been used since purchase or for more than 3 months.	Charge for at least 8 hours.
The water pressure is too low.	The tip is clogged.	Replace the tip.
	The filter is blocked.	Clean the filter.
The appliance can only operate for a short time after fully charged.	The battery has reached the end of its life.	Contact Service Hotline.
	Insufficient charging time.	Charge for at least 8 hours.
The charging indicator does not light on when charging.	The USB cable is not connected correctly.	Connect the power supply correctly.
	The charging circuit has a fault.	Contact Service Hotline.

### ONE YEAR LIMITED WARRANTY

This Conair Appliance is guaranteed for consumer use for one year. To obtain this warranty service, please keep the original purchase receipt. This guarantee is no longer valid in case of alteration or repair by any unauthorized persons and only covers defects under normal use from the date of purchase. If your product does not perform satisfactorily because of defects in materials or manufacture, it will be repaired or replaced under guarantee through the Conair Service Centre. Conair reserves the right to amend terms and specifications without prior notice.

### OUT OF GUARANTEE

Call Conair Service Centre in your region as listed in this Instruction Booklet for details of repair service. Or email to us at email\_asia@conair.com with your name, contact telephone number and email address and details of the fault.

### CONAIR CUSTOMER SERVICE:

35/F., Standard Chartered Tower, Millennium City,  
388 Kwun Tong Road, Kwun Tong, Kowloon, Hong Kong  
Conair Customer Hotline: (852) 2751 4604  
E-mail: email\_asia@conair.com  
Website: www.hkconair.com