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CWJ201H

IB-21/360

CONTINENTAL CONAIR LIMITED CHINA CREATIVE DEPARTMENT

Date	1/20/2022
Customer	VS-HK
Model No.	CWJ201H
Description	IB
Job No.	
Reference No.	IB-21/360
Revision No.	0
Color	1C+1C(Black)
Open Size	552mm(W) X 214mm(H)
Folded Size	92mm(W) X 214mm(H)
Material	80gsm Artpaper
Prepared by	Betty Bai
Approved by	Phoebe Laung/ Andy Huang
Designer	
Vendor	

折疊方法
Folding method:
12 pages (6 x 2)

CONAIR PORTABLE BATTERY DRIVEN WATER JET

INTRODUCTION

This portable water jet can improve your oral hygiene. It will effectively flush out the debris, reduce plaque between teeth, strengthen the gum and prevent tartar build-up.

SAFETY INSTRUCTION: IMPORTANT SAFEGUARDS

When using an electrical appliance, especially when children are present, the following safety precautions should always be taken: Polythene bags over product or package may be dangerous. To avoid danger of suffocation, keep the wrapper away from babies and children. The bag is not a toy.

READ ALL INSTRUCTIONS BEFORE USING KEEP AWAY FROM WATER

SAFETY PRECAUTIONS

This product is washable under water tap (secure battery cover before cleaning).

The above mark indicates that this product is suitable for use in a bath or shower.

- This appliance is not intended for use by children under 5.
- Do not use any attachments other than those provided with this product.
- Use this appliance only for its intended use as described within this instruction booklet.
- Unless otherwise specified by a dental professional, fill the water tank with lukewarm water only. Do not overfill the water tank. It is recommended for the temperature of the water not to exceed 40°C to avoid scalding.
- Do not direct the water jet under the tongue, into the ear, or other delicate areas. The product can produce pressures that may cause damage in these delicate areas.
- Do not use the appliance without a water jet tip.
- Some mouthwashes may corrupt this product. Please do not add mouthwash to this product.
- This appliance is a personal care product. It is not suitable for commercial use or shared use by multiple patients.
- Do not drop or insert any object into an opening.
- Do not place on any surface while it is operating.
- Do not leave unattended while it is operating.
- In the event of damage, discontinue use immediately and contact your dealer for instructions on returning it for examination or exchange.
- This appliance must be powered by safety extra-low voltage.
- Before cleaning, always unplug the appliance and secure the charging receptacle lid firmly.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure they do not play with the appliance.
- Do not use dishwasher or microwave oven to clean the jet tips and handle.


MEDICAL WARNING

- If you have any medical concerns, consult your doctor before using this product.
- Patients with cardiac disease should consult their doctor before using this water jet.
- Contact your dentist prior to use if you have periodontal disease or have had oral or periodontal surgery or major orthodontic treatment in the previous 12 months.
- When you start any new oral hygiene program, your gums may initially bleed slightly. Contact your dentist if excessive bleeding occurs or if bleeding continues after two weeks of use.
- Do not use mouthwashes or oral rinses that contains iodine, decolorizer or tea tree oil.
- Remove any oral jewelry prior to use of this product.
- Do not use if you have an open wound on your tongue or in your mouth.
- Make sure only the tip of the water jet comes contact to the teeth, gum or tongue. If any other part of the handle contacts them, it may cause damage to these areas. Be careful not to press the tip against the gum, especially dental patients.
- This product is designed for oral cleaning only. Do not use the unit for any other purpose.

IMPORTANT

- We recommend using longer-lasting alkaline batteries.
- If you are not going to use your water jet for a long time or if the batteries are dead, remove the batteries from the unit.
- Clean the battery contacts and those of the device prior to battery installation.
- Remove used batteries promptly.
- Always purchase the correct size and grade of battery most suitable for the intended use.
- Replace all batteries of a set at the same time.
- Do not mix old and new batteries. Do not mix alkaline, standard (carbon-zinc), or rechargeable (nickel cadmium or nickel metal hydride) batteries.


PRODUCT DESCRIPTION

OPERATING INSTRUCTIONS

BATTERY IN STALLMENT

- Check and ensure both the product and your hands are dry. Turn battery door lock as indicated in the figure below to open the battery door.

- Flip the battery door down as indicated in the figure below. Align the two AA batteries with the positive and negative positions indicated inside the battery door. Turn the battery door to close it and check whether it is firmly locked in place.

NOTE: Before each use, make sure the battery door is firmly closed and securely locked. If it is not securely locked, water may get into the product, which will cause the water jet to fail to work properly.

BEFORE USE

- Adjust the tip, then push the tip downwards until it is firmly attached to the handle (you will hear a click) (see the figure below). If you need to change the tip in the course of using, you should turn it off first. Simply press the changing button to change the tip.

TIPS: For best results, please replace the tip every 6 months. If the tip becomes loose or not "snapping into" the handle, you also need to replace the tip.

- Remove the water tank by sliding it down, hold the water tank vertically and fill with water. Then replace the water tank from the vertical direction onto the main body (you will hear a click).

NOTE: Some mouthwashes may corrupt this product, please do not add mouthwash to this product.


- Slide the switch up to choose mode:
Normal: suitable for normal cleaning of teeth.
Soft: suitable for users with sensitive teeth.

TIPS: In order to get a comfortable experience, when using the product for the first time, it is recommended to use soft mode for one week and change to normal mode until you get

adapted to this product.

USING THE PRODUCT



- Ensure the water jet is in off position.
- Hold the water jet vertically, place the tip into your mouth with the tip directing at spaces between your teeth or the gum. Keep your mouth slightly close to avoid water spraying outside and water easy running out of mouth.

- Slide the ON/OFF switch up to turn the water jet on.

NOTE:

- Direct the water jet at a right angle to the teeth or gum.
- Move along the teeth, clean the teeth in sequence.
- Direct the tip to the gumline at a right angle to the gum.
- Clean the teeth, braces, crowns and dentures thoroughly.
- Do not direct the water jet into tooth pocket.

- To turn the water jet off, simply slide the ON/OFF switch down.
- For the sake of your health, please follow the following steps after each use:
 - Open the water tank lid and pour out the water from water tank.
 - Press the power button. The water jet will start to work and drain out the water in the appliance.
 - Press the power button again. The flusher will shut off.
 - Replace the water tank lid.
 - Wipe off any water on the flusher with a clean and dry cloth.

DAILY CLEANING

NOTE:

- Before cleaning, please make sure that the battery cover is securely locked. If the battery cover is not locked firmly, the liquid may enter the water jet, and the water jet may be unable to charge or work after cleaning.
- Use clean water or neutral detergent to clean this product. Do not use corrosive or abrasive cleansers because they would cause damage to the product. Do not clean this product with hot water over 50°C.
- Do not use dishwasher or microwave oven to clean the tips, handle and USB cable.

- If the residues in the tip or water tank are not removed properly, there may be a hygiene problem.

CLEANING MAIN BODY:

- Clean the main body with a dry or damp cloth.
- Never immerse the main body into water.

CLEANING WATER TANK:

- Wipe the water tank with a dry or damp cloth, or remove it to rinse under running water.
- If you are not going to use it for an extended period, please pour out the water and let it dry naturally.

CLEANING WATER JET TIPS:

- Wipe with a dry or damp cloth, or remove it to rinse under running water. Then let it dry naturally.
- Do not bend, pull or twist the water hose.
- For best results, please replace the tip every 6 months. If the tip becomes loose or not "snapping into" the handle, you also need to replace the tip.

STORAGE

If you are not going to use the appliance for an extended period, in order to avoid bacteria breeding, please follow the following steps:

- Open the water tank lid, pour out the water.
- Press the power button, the water jet will start to work and drain out the water in the appliance.
- Press the power button again to turn off the flusher.
- Remove the water tank by sliding it down.
- Clean the handle, water tank, tips with a clean dry or damp cloth.
- Leave the handle, water tank, tips and charging base dry out of direct sunlight.
- Close the water tank lid and slide the water tank back to handle.
- Store the handle, tips and USB charging cable in a cool dry place.

TROUBLESHOOTING

Problems	Possible Reasons	Solutions
The water jet does not work	The water jet has not been used since purchase or for more than 3 months	Check whether batteries are installed in the appliance
The water pressure is too low	The tip is clogged	Replace the tip
	The filter is blocked	Clean the filter

DISPOSAL

This product requires dry batteries. When the batteries are dead or the product reaches the end of its useful life, please recycle or dispose the batteries according to your local waste regulations. If you have any questions, please consult the waste disposal office.

- In order to protect the environment, when the batteries run out, discard them in a recycling bin. If there is a leakage, be careful to avoid skin contact and discard the batteries in a recycling bin.
- Before discarding the product, ensure the batteries are removed and recycled or disposed according to your local waste regulations.
- Never discard the batteries into fire unless they are incinerated under controlled conditions. If this precaution is not followed, explosions may occur. Do not attempt to charge disposable batteries, because this may cause leakage, fire or explosion.

ONE YEAR LIMITED WARRANTY

This Conair Appliance is guaranteed for consumer use for one year. To obtain this warranty service, please keep the original purchase receipt. This guarantee is no longer valid in case of alteration or repair by any unauthorized persons and only covers defects under normal use from the date of purchase. If your product does not perform satisfactorily because of defects in materials or manufacture, it will be repaired or replaced under guarantee through the Conair Service Centre. Conair reserves the right to amend terms and specifications without prior notice.

OUT OF GUARANTEE

Call Conair Service Centre in your region as listed in this Instruction Booklet for details of repair service. Or email to us at email_asia@conair.com with your name, contact telephone number and email address and details of the fault.

CONAIR CUSTOMER SERVICE:

35/F., Standard Chartered Tower, Millennium City,
388 Kwun Tong Road, Kwun Tong, Kowloon, Hong Kong
Conair Customer Hotline: (852) 2751 4604
E-mail: email_asia@conair.com
Website: www.hkconair.com